

Compatibility of Bagnoli™ and Trigno™ Systems with External Hardware and Software Platforms

ACQUISITION SOFTWARE

- ✓ **EMMA**
- ✓ **LabChart**
- ✓ LabView
- ✓ MATLAB

BIO SIGNAL - DAQ HARDWARE

- ✓ ADInstruments PowerLab
- ✓ Biopac
- ✓ National Instruments

FORCE PLATES - TREADMILLS

- ✓ ADAL3D
- ✓ AMTI
- ✓ Bertec
- ✓ Kistler
- ✓ Motekforce Link
- ✓ Treadmetrix

ISOKINETIC MACHINES

- ✓ Biodex
- ✓ HUMAC

PRESSURE MEASUREMENT SYSTEMS

- ✓ Novel
- ✓ RsScan
- ✓ Tekscan

MOTION CAPTURE

- ✓ C-Motion
- ✓ Codamotion
- ✓ Contemplas
- ✓ Dartfish
- ✓ Motion Analysis Corporation
- ✓ **The MotionMonitor**
- ✓ NDI Optotrak
- ✓ Organic Motion
- ✓ Optitrak
- ✓ Phasespace
- ✓ Qualisys
- ✓ Simi
- ✓ STT Systems
- ✓ TEA
- ✓ Vicon
- ✓ Xsens

OFF-THE-SHELF CAMERA


- ✓ Basler


NEUROPHYSIOLOGY DEVICES


- ✓ Ant Neuro
- ✓ CED
- ✓ G.Nautilus
- ✓ Magstim

OTHER INSTRUMENTS

- ✓ Biometrics
- ✓ Interface Force
- ✓ Interlink Electronics


 Acquisition Software > **LabChart**





 Acquisition Software > **EMMA**

 Motion Capture > **The MotionMonitor**


 Compatibility refers to software integration and/or system synchronization. Additional devices and software may be supported. For full details, please contact support@delsys.com.



Compatibility of Bagnoli™ and Trigno™ Systems with External Hardware and Software Platforms

For more details on specific integration information, please contact delsys@delsys.com.

	BAGNOLI™ DESKTOP EMG SYSTEM	TRIGNO™ WIRELESS SYSTEM (ANALOG + DIGITAL)	TRIGNO™ WIRELESS SYSTEM (DIGITAL)
TRIGGERING	Connect Trigno™/Bagnoli™ + 3 rd party hardware to Delsys Trigger Module for data synchronization.		
ANALOG	BNC or 68-pin analog outputs enable direct connections into variety of acquisition hardware	Optional 68-pin analog outputs enable direct connections into variety of acquisition hardware	
DIGITAL		Output data digitally into specialized software via SDK Combine with triggering for complete synchronization	
OTHER		Trigno™ Auxiliary Sensors designed to collect specific 3 rd party signals directly through Trigno™ system. Compatible with load cells, goniometers, footswitch pressure leads, dynamometers, custom analog devices. <div style="display: flex; justify-content: space-around; margin-top: 20px;"> <div style="text-align: center;">  <p>Load cells</p> </div> <div style="text-align: center;">  <p>Footswitch pressure leads</p> </div> </div>	



23 STRATHMORE ROAD
NATICK, MA 01760 • USA

DIRECT 508.545.8200
FAX 508.975.4551

delsys@delsys.com
WWW.DELSYS.COM